
Invia® Y-Connector with Quick-Connector for Negative Pressure Wound Therapy

EN INSTRUCTIONS FOR USE Invia® Y-Connector with Quick-Connector

Intended use

The Invia® Y-connector with Quick-connector is intended for use with the Invia Negative Pressure Wound Therapy (NPWT) system to connect two devices (Invia® FitPad and/or Invia Drain adapter) to one Invia pump.

The Invia® Y-connector is designed for use on a wide range of wounds indicated for NPWT. For indications, contraindications, warnings and cautions of the Invia® Liberty™ and Invia® Motion™ NPWT systems consult the Invia Liberty instructions for use or Invia Motion instructions for use.

To consider before use

The Invia NPWT system is intended for use in hospitals, clinics, Long Term Care (LTC) and Home Care (HC) settings by qualified and trained medical personnel adhering to the instructions for use. To apply, use clean/aseptic or sterile techniques in accordance with local protocol.

Application instructions

1. Connect the Invia FitPad and/or the Invia Drain adapter to the Y-connector.
2. Connect the Y-connector to the canister tubing and turn on the pump.

Recommended pressure level is the same as when dressing a single wound.

Frequency of change and disposal

The Y-connector can be considered as an extension of the canister tubing. The Y-connector should be changed when the canister tubing is changed, but no less than once a week. Dispose the Y-connector in accordance with local environmental regulations.

Monitoring of the tubing

When two Invia FitPads are used, only the FitPad tubing connected to the straight arm of the Y-connector (A) is flushed with air. A blockage in the FitPad tubing connected to the curved arm of the Y-connector (B) will not be detected by the pump.

1. Instruct the patient to frequently check tubing connected to the curved arm of the Y-connector (B) to make sure tubing is not kinked or blocked.
2. If blockage is identified along FitPad tubing connected to the curved arm of the Y-connector (B), patient is instructed to switch the Invia NPWT pump off and contact the healthcare professional.

Warnings

- When two Invia FitPads are used, only the FitPad tubing connected to the straight arm of the Y-connector (A) is flushed with air. A blockage in the FitPad connected to the curved arm (B) will not be detected by the pump.

1. Instruct the patient to frequently check tubing connected to the curved arm of the Y-connector (B) to make sure tubing is not kinked or blocked.
2. If blockage is identified along FitPad tubing connected to the curved arm of the Y-connector (B), patient is instructed to switch the Invia NPWT pump off and contact the healthcare professional.

- Do not connect wounds with different etiology in which cross contamination may occur.

- Avoid using a Y-connector to connect wounds that would be optimally treated with different pressure settings.

- It is not recommended to use a Y-connector to connect grafts and/or flaps.

- Do not reuse the Y-connector. If reused performance of the product may deteriorate, cross contamination may occur.

- The product is sterile. Do not use if inner package is damaged or opened prior to use. Do not re-sterilize.

Cautions

- Only one Y-connector should be used for each pump unit.

- If there is a risk for imprints from the Invia Y-connector, place a skin protector between skin and Y-connector.

- Invia Y-connector with Quick-connector is only compatible with Invia NPWT systems with Quick-connector.

CAUTION: Federal US law restricts this device to sale by or on the order of a physician.

	This symbol indicates the device is sterilized using ethylene oxide.		This symbol indicates to not use the device if package is damaged.
--	--	--	--

	This symbol indicates a prescription device. CAUTION: Federal US law restricts this device to sale by or on the order of a physician (for US only).		This symbol indicates the name and the address of the manufacturer.
--	---	--	---

	This symbol indicates a single use device. Do not reuse the device.		This symbol indicates the date of manufacture.
--	---	--	--

	This symbol indicates manufacturer's catalog number.		This symbol indicates manufacturer's batch code.
--	--	--	--

	This symbol indicates to follow the instructions for use.		This symbol indicates the temperature limitation for operation, transport and storage.
--	---	--	--

	This symbol indicates that the device should not be used after the date shown.		This symbol indicates the atmospheric pressure limitation for operation and storage.
--	--	--	--

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