Invia[®] Liberty[™]

NEGATIVE PRESSURE WOUND THERAPY SYSTEM

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Assistance

Read the entire instruction manual before trying to operate the device.

If you have any questions about your Invia Liberty Negative Pressure Wound Therapy System, please contact your healthcare provider/physician at the number below:

HEALTHCARE PROVIDER/PHYSICIAN
CONTACT INFORMATION:

Keep this Patient Instructions for use booklet in an easily accessible location.

Table of Contents

EN Instructions for Use	3
ES Instrucciones para uso	29
FR Mode d'emploi	58
General Warnings.....	4
Purpose of the Device.....	7
Intended Use	7
When Not to Use Device (Contraindications).....	8
Description of the Device.....	9
Display.....	10
Conditions that Affect Use	10
General Safety at Home.....	11
Sleeping.....	12
Showering, Bathing and Washing.....	13
When You Leave the House	13
Set Up Instructions	14
Set Up Carrying Case.....	16
Charge Battery.....	17
Change Canister.....	19
Notifications and Troubleshooting.....	21
Alarm table.....	22
Cleaning/Disinfection.....	26
Guarantee.....	26
Maintenance/Safety related Check.....	26
Disposal.....	26
Technical Specifications.....	27
Signs and Symbols.....	28
Electromagnetic Emissions.....	86

General Warnings

WARNING

The individual responsible for operating this device on a daily basis (patient, caregiver) must be able to read, understand and follow directions provided in this IFU, on the pump display and by the healthcare professional. The individual must be able to hear and respond to auditory signals.

Only use the Invia Liberty NPWT System on the person for whom it was ordered and only for its intended use.

Do not perform the therapy without first consulting with healthcare personnel and caregivers trained on the instructions for use and the healthcare professional's supervision.

No modification of this equipment is allowed.

WARNING

Visually check at the wound dressing and Invia Liberty frequently.

If you see large amounts of blood in the dressing, tubing or canister:

- 1 Turn off the pump by pressing  and hold for 3 seconds.
- 2 Apply pressure to your wound.
- 3 Seek immediate Emergency Medical assistance.

If you have difficulty breathing, seek immediate Emergency Medical Assistance.

 **WARNING**

Immediately contact your healthcare professional if you experience the following symptoms:

- Significant vision or hearing loss
- Pain, sore or red area around the wound
- Swelling, itching, or rash around the wound
- Fever
- Strong odor from the wound
- Thick green drainage or any change in color of drainage
- Nausea or vomiting or diarrhea
- Sore throat with swelling
- Dizziness or confusion
- Allergic reactions like redness, rash or hives, swelling or severe itching

 **WARNING**

Once the healthcare professional has set the pump pressure and mode settings do not press any other button than the on/off button [], the standby button [] or the mute button [].

 **WARNING**

The pump must be used 24 hours per day. Do not stop the pump unless instructed by your healthcare professional or if excessive bleeding occurs.

If the pump is stopped for more than the time frames shown below, the dressing must be changed. Contact your healthcare professional to request a dressing change.

Gauze dressings: Change dressing if the pump is stopped more than 24 hours.

Foam dressings: Change dressing if the pump is stopped more than 2 hours.

 **WARNING**

If the pump stops working and cannot be switched on, immediately call your healthcare professional.

Purpose of the Device

Your healthcare provider has prescribed Negative Pressure Wound Therapy (NPWT) for you.

Negative Pressure Wound Therapy has been shown to help promote healing in several different kinds of wounds. Negative pressure (suction) is delivered to a wound through a pump. The suction, or negative pressure, helps to remove fluid and close your wound. Your healthcare professional will place a special dressing on your wound and a tube will go from your wound to the canister on the pump. The pump is then turned on to the prescribed pressure to remove the fluid from your wound. The fluid then goes into the canister. The treatment will be applied 24 hours a day.

Intended Use

The Medela® INVIA Liberty Wound Therapy is indicated to help promote wound healing, through means including drainage and removal of infectious material or other fluids, under the influence of continuous and/or intermittent negative pressures, particularly for patients with chronic, acute, traumatic, subacute and dehisced wounds, partial-thickness burns, ulcers (such as diabetic or pressure), flaps and grafts.

When Not to Use Device (Contraindications)

Contraindicated for patients with:

- Malignancy of the wound
- Untreated osteomyelitis
- Unmanaged malnutrition
- Non-enteric fistula
- Unexplored fistula
- Necrotic tissue with eschar present
- Do not place Invia Wound Therapy dressing over exposed blood vessels or organs
- Do not place directly over anastomoses or sutured vessels

If you have concerns, speak with your healthcare professional.


Description of the Device

Check the delivery package of the Invia Liberty NPWT system for completeness and general condition.


WARNING

Do not use any other accessories than the ones described below.

Display

Buttons to Operate the Invia Liberty

Conditions that Affect Use

CAUTION

The Invia Liberty must remain in an upright position during use.

General Safety at Home

2

Please read and observe these warnings and safety instructions before operation. These Instructions for Use must be kept with the device. Please note that these Instructions for Use are a general guide for the use of the product. Medical situations must be addressed by a physician.

CAUTIONS/WARNINGS

- Do not use extension cords with the Invia Liberty.
- Position the power cord so that it will not cause someone to trip and fall.
- Never place the power supply cord around your neck.
- Keep power supply cord away from hot surfaces.
- Do not unplug the power supply cord from the wall outlet by pulling on the cord.
- Do not try to fix or repair the Invia Liberty. Contact your healthcare professional if the pump is not working.
- Keep the Invia Liberty clean and dry.
- Never place the Invia Liberty in water or liquids.
- If the Invia Liberty gets wet, dry with a towel.
- Do not dry the Invia Liberty in a microwave.
- Keep the Invia Liberty away from children and pets.
- Supervision is necessary when the Avance wound-pump is used in the vicinity of children.
- Keep the Invia Liberty protected from direct sunlight.

- The use of mobile telephones, local area networks (LAN) including wireless (WLAN), walkie-talkies (two-way radios) telephone sets can affect the Invia Liberty. This may cause insufficient battery charging or may lead to an internal fault. See the chapter Notifications and Trouble-shooting for actions to take in the case of an internal fault. Certain mobile phones may interfere with the Invia Liberty NPWT system when receiving and making calls close to the pump. Prior to using a mobile phone near the Invia Liberty NPWT system, check that it continues to function while making or receiving a call.

If you observe a change in the vacuum level or an error condition, move the Invia Liberty NPWT system at least 1 meter (3.3 feet) from the mobile phone when making or receiving calls. Turn the Invia Liberty NPWT system on and off again to reset it after an error condition. Whenever possible, a safety distance of min. 3.3 feet (1 meter) between the Invia Liberty and a mobile phone is recommended.

Sleeping

- Make sure the Invia Liberty tubing will not become kinked.
- Prevent the Invia Liberty from falling to the floor during sleep.
- Always charge the Invia Liberty when going to bed.

Showering, Bathing and Washing

23

- Do not use the Invia Liberty while bathing or showering.
- The Invia Liberty Pump can be disconnected from the dressing for a short period of time (up to 20 minutes) in order to change clothes or perform personal hygiene. Discuss bathing plans with your healthcare professional.
- Prevent the Invia Liberty from getting wet.
- Do not put the Invia Liberty where it can fall into a shower, sink or tub.
- Do not touch the Invia Liberty when it has fallen into water. Unplug device immediately from electrical outlet and disconnect the pump from the dressing.

When You Leave the House

- Check the battery status and recharge if necessary before leaving the house.
- For your convenience, place the Invia Liberty into the carrying case.
- Don't forget to take a spare Canister and power supply with you.

Set Up Instructions



WARNING

Check the Invia Liberty frequently to make sure it is working. Visually check the wound dressing, area around the wound and fluid in the canister for signs of infection or other complications as provided by your healthcare professional.

Turn ON

Switch Invia Liberty Pump ON by pressing [].

Display view when the pump is running. Run/Standby indicator is rotating [].

WARNING

Do not try to perform therapy without the supervision of your healthcare professional.

**WARNING**

NPWT must run 24 hours a day. Do not stop the therapy unless instructed by your healthcare professional or if excessive bleeding occurs.


**WARNING**

The dressing should be changed and therapy restarted by the healthcare professional if therapy is discontinued for:

- More than two (2) hours when using foam.
- More than 24 hours when using gauze.

Turn OFF

Switch Invia Liberty Pump OFF by pressing [] for 3 seconds. Press only once.

Set Up Carrying Case

The Invia Liberty is intended to be worn in a carrying case during operation.


- 1 Open the carrying case.
- 2 Slide the Invia Liberty into the designated pocket.
- 3 Close the velcro of the pump.

The shoulder strap can also be used as a belt strap. To accomplish this, release the strap at the sides of the carrying case and pull it through the two loops located on the back of the case.


Charge Battery

When the battery is fully charged it should last approximately 10 hours. If the Invia Liberty is running continuously due to an air leakage, the battery may need to be recharged earlier (after approximately 4 hours).

Battery Status Symbols

If the battery becomes empty, please follow the instructions below.

Open the electrical outlet port. Connect the Invia Liberty to the supplied electrical plug from the power supply provided with the pump. Plug the other end of the power supply into a wall outlet. Moving bars on the battery display will indicate charging is occurring.

The battery will be fully charged after approximately 4 hours. Once charged  will appear constantly on the display. After charging is completed, unplug the power supply from the wall and disconnect the Invia Liberty from the power supply.

Close the electrical outlet port by pressing the protection cover down.

CAUTION

NPWT therapy must be continued while charging the pump's battery. Do not turn the pump off while charging.

Always charge the battery over night.

Change Canister

Change Canister when it appears full or when the alarm sounds, in accordance with instructions on the pump display.

1. Wash hands.
2. Clamp double lumen tubing.
3. Press the Standby button [ > 3 seconds] and the pump motor will stop running.

4.  Press release button at top of pump and remove canister.

5.  Seal used canister with cap.

6.

Unpack new canister, position the feet at the bottom of the canister into the grooves at the base of the pump. Push the canister up and click into the Invia Liberty Pump.

CAUTION

Do not use an Invia Liberty canister if the sterile packaging is damaged.


7.

Switch the Invia Liberty on by pressing [] and slide the pump into the carrying case.


8. Discard the used canister according to instructions provided by your healthcare professional.

Notifications and Troubleshooting









Don't panic if an acoustic signal sounds. Follow the instructions below.

- 1 If the Invia Liberty is in the carrying case open the lid of your carrying case.
- 2 Press [] to mute the acoustic signal for 60 seconds (can be pressed every time the signal sounds, with the exception of "Battery Empty").

Alarm table

	Fault number	Problem description on the display	Troubleshooting on the display	D - - C - C - -
Alarm	301			D - - C - C - -
	305			I (I
	306			C T
	313			C

Remarks/potential cause of fault

Dressing:

- Check dressing for air leakage. Press firmly around the edges of the dressing, around the drain tube or on the Transfer Pad.
- Apply some additional film dressing to seal the leaking area.

Connectors:

- Ensure that the tube connected to the dressing is connected properly to the Canister tube.
- Ensure that the canister tube is inserted straight into the pump.

Canister:

- Ensure that the canister is properly inserted, release the canister and reposition.
- Ensure that the O-ring/gasket, placed beside the canister tubing on the pump is not missing. Additional O-ring is available via Medela representative.

Recharge the battery either by placing the Invia Liberty Pump in the Docking Station or plug in the charger to the electrical outlet port on the pump.

Remaining time of battery is approximately 15 minutes.

Change the canister, see chapter “Change Invia Liberty Canister and Invia Liberty Tubing”.

Change canister, see chapter “Change Invia Liberty Canister and Invia Liberty Tubing”

	Fault number	Problem description on the display	Troubleshooting on the display
Warning	401		
	405		
Internal fault			

CAUTION

If you are unable to identify any symbol or don't understand the alarm on the display, switch off the pump by pressing [] for 3 seconds. Contact your healthcare professional.

Remarks/potential cause of fault
<p>Recharge the battery either by placing the Invia Liberty Pump in the Docking Station or plug in the charger to the electrical outlet port on the pump.</p> <p>Remaining time of battery is approximately 30 minutes.</p>
<p>If the pump is in Standby mode for more than 5 minutes, an alarm will go off.</p> <p>To continue therapy press “On” [] or switch off the pump by pressing [ > 3 seconds].</p>
<p>Restart the pump. If internal fault remains, turn off by pressing [ > 3 seconds] and contact Medela Customer Service.</p>


WARNING

If the pump fails or cannot be switched on call your healthcare professional immediately.

Cleaning / Disinfection

The Invia Liberty can be wiped with a damp cloth using a mild household cleaner.

Guarantee

Warranty period is for 2 years after date of delivery. The manufacturer is not liable for any damage or consequential damage caused by incorrect operation, inappropriate usage as well as use by unauthorized persons.

Maintenance / Safety related Check

Service work may only be carried out by authorized personnel. A safety-related check are confined to visual inspection of the housing and charger for damage and must be performed prior to each use.

If Invia Liberty Pump has not been in use, the battery must be charged approximately once every 6 months to ensure optimum function.

Disposal

Ask your healthcare professional how to dispose of used Canister at the end of your therapy.

Technical Specifications

2

Accessories Storage Conditions

Accessories must remain in the packaging and stored at a temperature range of -4° Fahrenheit and +122° Fahrenheit (-20° Celsius and +50° Celsius).

Operation Environmental Conditions

The Invia Liberty and accessories must be operated within a temperature range of +41° Fahrenheit and +104° Fahrenheit (+5° Celsius and +40° Celsius). Do not operate the products in extreme cold or heat. Avoid high relative humidity (15–93% is permissible).

CAUTION

Do not operate the Invia Liberty above 9,842 feet above sea level (3000 meter above sea level).

Protection Class IP33

The Invia Liberty is protected against ingress of solids larger than 0.1 inch (2.5 millimeter) and dripping water.

Signs and Symbols

This Symbol indicates a CAUTION or WARNING associated with the Device (see page 3).

This Symbol indicates a class II Device.

Do not dispose the device together with unsorted municipal waste (for EU only).

This Symbol indicates the Date of Manufacture (four Digits for the Year and two Digits for the Month).

This Symbol indicates the Name and the Address of the Manufacturer.

This Symbol indicates the Device is sterilized using ethylene oxide.

This Symbol indicates that this is a Prescription Device. Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner. (for US only).

This Symbol indicates a type BF applied Part.

This Symbol indicates manufacturer's Catalog Number.

This Symbol indicates manufacturer's Serial Number.

This Symbol indicates manufacturer's Batch Code.

This Symbol indicates the protection against the ingress of solid foreign objects and against harmful effects due to the ingress of water.

This Symbol indicates that the Device should not be used after the End of the Year and Month shown.

This Symbol indicates to follow the Instruction for use.

This Symbol indicates the Humidity Limitation for Operation, Transport and Storage.

This symbol indicates the compliance with additional USA and Canada safety requirements for medical electrical equipment.

This Symbol indicates the atmospheric Pressure Limitation for Operation, Transport and Storage.

This Symbol indicates to not use the Device if package is damaged.

This Symbol indicates to keep the Device dry.

This Symbol indicates the number of items n that the content is sufficient for.

This Symbol indicates to keep the Device away from sunlight.

pcs

This Symbol indicates the number of items.

This Symbol indicates to handle the fragile Device with care.

This Symbol indicates the Direct Current Socket.

This symbol indicates the compliance with the essential requirements of the Council Directive 93/42/ EEC of 14 June 1993 concerning medical devices.

This symbol indicates a Single Use Device. Do not reuse the Device.

This Symbol indicates the Temperature Limitation for Operation, Transport and Storage.

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